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[www.proactive-painmanagement.com](http://www.proactive-painmanagement.com)

## INSTRUCTIONS: BACK WRAP

Back wrap should be applied for 20 minutes, or as recommended by a doctor. Remove and wait at least 20 minutes. Repeat as needed. To further reduce hot or cold level, place wrap in soft cloth covering before applying.

### For Heat Therapy:

Microwave wrap per the below suggested heating times and test wrap temperature. Heat for additional 5 seconds as needed until desired temperature is reached.

#### Suggested heating time for one (1) back wrap:

700 Watts 65 seconds	1000 Watts 55 seconds	1250 Watts 45 seconds
-------------------------	--------------------------	--------------------------

#### Do not exceed heating time for one (1) back wrap:

700 Watts 85 seconds max.	1000 watts 75 seconds max.	1250 Watts 65 seconds max.
------------------------------	-------------------------------	-------------------------------

### For Cold Therapy:

Freeze for at least 2 hours in a freezer storage bag.

### Caution:

Always test back wrap temperature prior to use. Overheating can cause damage or destroy it. Prolonged use may cause injury to user including burns. For external use only. Do not use on infants. Never sit, lie, or lean against the wrap. Do not puncture wrap or ingest contents. Do not apply to sensitive/insensitive skin or open wounds. Use only for intended purposes. Do not fall asleep using wrap. Use only under adult supervision. For chronic, extreme or prolonged pain, please contact your physician.

### Storage:

Always keep wrap in a freezer storage bag when kept in the freezer.

### Cleaning:

To clean wrap, use mild soap and water and air dry before reusing.

Manufactured and printed in China for:

A.M.G. Medical Inc.

Montréal, QC H4T 1V5 1-800-363-2381 West Chazy, NY 12992 1-888-412-4992

[www.proactive-painmanagement.com](http://www.proactive-painmanagement.com)

Rev 1 0417 P/N 911-740804