


Airgo[®]

Flex-Grip[™]

Cane Tip




Instructions:

1. Remove the old tip from your cane.
2. Slip on the Flex-Grip[™] cane tip, making sure that the AMG logo is aligned with the front of the cane (the side that is facing forward when you walk). This will ensure optimal grip and traction for a more natural walking motion.

Fits most aluminum canes

To verify that the Airgo[®] Flex-Grip[™] cane tip is right for your cane, remove the old cane tip and match the shaft end with the circle below.



Questions? Comments?

www.amgmedical.com • 1-800-363-2381



Model no. 735-001
Made in China for AMG Medical Inc.
Montréal, QC H4T 1V5

REV 1 0206
P/N 815-381