


### 700-983 Hugo® Sidekick™ Side-Folding Rolling Walker with a Seat Tangerine



**700-983**  
**Hugo® Sidekick™ Side-Folding Rolling Walker with a Seat**  
**Tangerine**

Item #	AMG Part #	Item Description	Qty total	Qty sold
1.	<b>C23-00016</b>	Handle Bars with Brake Handle Assy (Handle Bar, Brake Handle with Housing & Hardware, Left / Right Ortho-handle Grip & End Plug) - L & R	1 PR	1 PR
2.	<b>C23-00017</b>	Vinyl Seat 18" x 9" and Seat Strap with Hardware (2 steel bars & 4 Cap Screws)	1 ST	1 ST
3.	<b>C23-00018</b>	Brake Cable and Pad Assy (Cable, Sleeve, Brake Pad L&R, Cable Adjustment Mechanism, Spring with Spring Sleeve, Cable Stopper with Set Screw & Fishtail Stopper) – L & R	1 PR	1 PR
4.	<b>C23-00019</b>	Mounting Brackets with Hardware - L & R	1 PR	1 PR
5.	<b>C23-00020</b>	6" Rear Wheel (3 Spokes) with Hardware (Curved Washer & Locknut M8) – 2 ea / set	1 ST	1 ST
6.	<b>C23-00021</b>	Backrest for Sidekick	1 EA	1 EA
7.	<b>C23-00022</b>	Set of 2 E-Clips dia=22mm, 2 E-Clips Dia=19mm, 2 Tightening Knobs & 2 Hex Bolts M6x32mm, 2 Cable Clips, 4 Wheel Caps & 2 Fork Caps	1 ST	1 ST
8.	<b>C23-00024</b>	Front Leg (Tangerine) with Fork, Wheel and Hardware – L & R	1 PR	1 PR
9.	<b>C23-00025</b>	Hugo-Sidekick under Seat Soft Bag, Black	1 EA	1 EA