


770-302 CLINELL COMMODE



Item #	AMG Part #	Item Description	Qty Per Unit	Pkg Qty
1.	A15-411	Backrest, White Plastic	1 EA	1 EA
2.	B15-401	Seat with Lid, White Plastic	1 ST	1 ST
3.	B15-351	Footrest Left & Right, Grey Plastic	1 PR	1 PR
4.	B15-101	Caster with Locking System, Grey Plastic (4 ea / set)	1 ST	1 ST
5.	B15-251	Armrest Left & Right, White Plastic	1 PR	1 PR
6.	B15-501	Reusable and Disposable Bedpan Support, White Plastic	1 ST	1 ST