

**NEW!**

# Airgo™

## Lightweight Wheelchair

*For who you are, where you live, what you need...*



The Airgo™ by AMG Medical is a high performance, lightweight wheelchair that combines versatility, simplicity, durability along with superior value. Available in a range of sizes and a selection of options to fit the needs of most users.





# Lightweight Wheelchair

*For who you are, where you live, what you need...*

## Features

- High strength powder coated carbon steel frame for durability and low maintenance.
- Dual axle positions allow for varying seat-to-floor heights.
- Urethane rear tires and front casters for lighter weight and enhanced performance.
- Maintenance free dual-sealed precision fork, caster and wheel bearings.
- Dual armrest-cradle design for both flip-back and easy removal.
- Heavy duty, cushioned seat and back with anti-sag construction for added comfort and durability.
- Flame-retardant nylon upholstery resists cracking and UV damage.
- Seat and back screw inserts for longer life, even after multiple upholstery changes.
- Frame design accommodates custom backs and seating.
- Meets RESNA standards for manual wheelchairs in product performance.

## Available Accessories

- 701-300 Airgo™ Back Pack
- 701-310 Rear Anti-tip Device
- 701-320 Brake Handle Extension
- 701-340 Oxygen Tank Holder
- 701-350 Telescoping IV Pole

*See just how easy it can be to configure your Airgo™ at [www.amgmedical.com/airgo](http://www.amgmedical.com/airgo)*

## Specifications

<b>SEAT WIDTHS:</b>	16", 18", 20", 22"
<b>SEAT DEPTHS:</b>	16", 18"
<b>SEAT-TO-FLOOR HEIGHTS:</b>	15.5", 17.5", 19.5"
<b>BACK HEIGHTS:</b>	16"-17"-18" adjustable

### FRONT RIGGINGS:

- Swing-away footrests
- Swing-away padded elevating legrests

### ARM STYLES:

- Full-length, fixed height flip-back, removable
- Desk-length, fixed height flip-back, removable
- Full-length, adjustable height flip-back, removable
- Desk-length, adjustable height flip-back, removable

**PRODUCT WEIGHT:** 15 kg (33 lbs)

**WEIGHT CAPACITY:** 250 lbs (16", 18" seat width)  
350 lbs (20", 22" seat width)

### WARRANTY:

- Lifetime on frame and cross-braces
- 6 months on upholstery and all other components



To learn more about Airgo™, contact your AMG Representative or call toll-free at 1-800-363-2381.



**AMG Medical Inc.**  
Montreal, QC H4T 1V5  
[www.amgmedical.com](http://www.amgmedical.com)

