

A drug-free way to control your pain!

ProActive™

Make Your Life Better



Electrotherapy Products

EASY • EFFECTIVE • DISCREET • SAFE



What is pain?

Pain is debilitating. Our quality of life is affected when we suffer from pain. It is our body's natural mechanism to warn us against potential injury or is a signal that other underlying conditions may exist.

There are 2 types of pain: acute and chronic.

Acute pain is pain that begins suddenly. It can be sharp or severe, and is the type of pain that can occur after surgery, dental work, strains, or sprains. Acute pain is short-term and disappears after treatment or when the body's natural healing powers resolves it.

Chronic pain is pain that continues after your body has healed from an illness or injury. There may be no identifiable cause for the persistent pain. It lasts for a longer period of time and interferes with daily life at all levels.

What is TENS?

TENS is a drug-free way to control acute and chronic pain. It is a form of electrical stimulation that is applied using electrodes placed upon the skin. The resulting action is to modulate the perception or feeling of pain.

TENS devices use two approaches to pain relief:

1. TENS devices target your sensory nerves, stimulating them to block pain signals and prevent their transmission to the brain.
2. TENS devices promote production of endorphins, the body's natural pain relieving substances.

What conditions would benefit from TENS?

- Arthritis
- Tendonitis/Bursitis
- Radiculopathy
- Phantom limb pain
- Spinal cord disorders
- Low back pain
- Carpal tunnel syndrome
- Peripheral nerve injuries
- Post operative pain
- Fibromyalgia

It is recommended to consult your health professional to determine if TENS is a suitable therapeutic device for your health needs.



ProActive™ TENS (715-500)



A drug-free way to control your pain!

- A CSA certified compact easy-to-use dual channel TENS nerve stimulator for active pain control that offers variable pulse width and frequency settings.
- Choice of 3 versatile stimulation modes: burst, normal and modulated waveform settings to meet most treatment protocols.
- Choice of 3 different timer settings: 30 minute, 60 minute and continuous for personal customization.
- 2 lead wires with built-in guard to protect against wire splitting.
- 4 easy to apply 2" round re-usable self-adhesive comfort gel electrodes
- On-the-go pouch pack
- Easy to understand instruction manual



ProActive™ Self-Adhesive Electrodes



100% Latex Free!

Item No.	Size	Shape	Qty/Pk
715-710	2"	Round	4
715-715	2"	Square	4
715-720	1.5" x 3"	Rectangular	4
715-725	4" x 6"	Butterfly	1

For use with most TENS and muscle electro stimulators

- Flexible cloth backing allows the electrode to easily and durably conform to all body surfaces
- Pre-wired and pre-gelled for easy and convenient application
- Long lasting, hypo-allergenic gel keeps skin soft
- Rounded edges for prevention of electrode curl-up
- Anchored pigtail prevents wire pull-out and increases usage life
- Available in a variety of shapes and sizes for application to multiple treatment zones

More TENS kits also available!



TENS Plus Dual Channel TENS kit with battery charger (715-400)

TENS Plus Dual Channel TENS kit with battery charger and re-usable self-adhesive gel electrodes (715-410)



A full line of accessories



Item No.	Description	Size	Qty/Pk
715-510	Carbon electrodes	2" Square	4
715-525	Self adhesive patches to attach carbon electrodes	3" Square	20
715-530	Lead wires for TENS and muscle electrostimulators		2
715-402	Rechargeable 9 Volt battery		1
715-405	9 Volt battery charger		1
715-550	Hypo-allergenic, water soluble electrode gel for use with carbon electrodes	60g Tube	1
715-706	Convenient on-the-go pouch pack for TENS devices		1

You can learn more about the ProActive line of electrotherapy products by calling 1-800-363-2381 or by visiting us at www.amgmedical.com

AMG AMG Medical Inc.
Montreal, QC H4T 1V5
www.amgmedical.com



REV. 1.0 0704
P/N 915-300