

# MedPro®

The Professional Choice



## Practical, hygienic and infection control friendly

- Less risk of cross contamination – single use vessels.
- Can be used for many applications.
- Soft surface, smooth and warm to the touch, for patient comfort.
- Will not adhere to moist skin.
- Easy to handle, easy to carry and light.
- Stackable – easy to store and separate.



**Bedpan Liner**  
762-640



**Lid for Bepan Liner**  
762-641



**Fracture Bedpan**  
762-653



**Lid for Fracture Bedpan**  
762-654



**Large Bowl (2L)**  
762-678



**Lid for Large Bowl**  
762-679



**Kidney Dish**  
762-650



**General Purpose Bowl (3 L)**  
762-700



**Small Bowl (1.3 L)**  
762-676



**Female Urinal**  
762-655



**Male Urinal**  
762-660



**Urine Specimen (1 L)**  
762-701

- Fluid resistance standard: PAS29:1999.
- Resistance / impermeability testing up to 24 hours (water: 24 hours, phenol: 20 hours, alcohol: 80 minutes.)
- Resistance temperature 35°C ± 3°C.

- Molded pulp is made from virgin pulp or recycled paper and cardboard waste, rosin, and water.
- Non sterile bulk packaging (recyclable cardboard case).

**Minimal impact on the environment:**

- Natural, biodegradable, recyclable, non-toxic (even when burnt), water proof, oil resistant, static free, shock absorbent, microwaveable.
- Suitable for disposable in most macerators, if facilities are available.

Rev1 0910

# MedPro®

The Professional Choice



## Disposable recycled paper containers

**Less time spent cleaning means more time for patient care.** With this in mind, MedPro® offers a wide range of disposable containers that make human waste disposal more cost-efficient and allow staff to maximize their time with patients. By reducing contact with human waste, they also help to reduce the risk of cross-contamination and help support your infection control program.



Containers made of 100% recycled paper.



**Institutional**















**Homecare**



**Private**

For more information about MedPro® products, contact your AMG representative or our customer service at 1-800-363-2381.

**AMG** AMG Medical Inc.  
Montreal, QC H4T 1V5  
www.amgmedical.com

Product	Description	Use With	Size	Maximum Capacity	N° of Units per Case	Case Size (L xW x H)	AMG N°	UPC
	<b>Bedpan Liner</b> Standard hospital bedpan liner. Designed for patient comfort and for use with either a bedpan support or commode. Will not adhere to moist skin. Easy to handle: flat bottom and end surfaces, large tick rim. Easy to carry.	MedPro Euro and AquaCare bedpan support seat (770-380)	14.1" x 12" x4" (35.9 x 30 x 10 cm)	2 L	100	18.1" x 12.2" x 14.75" (46 x 31 x 37.5 cm)	762-640	775757626406
	<b>Lid for Bedpan Liner</b> —	Bedpan Liner	11.5" x7" + tab (29 x 18 cm + tab)	—	200	19.3" x 12.8" x 7.5" (49.5 x 32.5 x 19 cm)	762-641	775757626413
	<b>Fracture Bedpan</b> Designed for patient comfort and for use with a slipper pan support. Slide under a less mobile patient. Will not adhere to moist skin. Easy to handle: flat bottom and end surfaces. Lightweight and easy to carry.	Slipper pan support	13.4" x 11.6" x 2.75" (34 x 29.5 x 7 cm)	1.3 L	100	13.6" x 12.4" x 14.6" (34.5 x 31.5 x 37 cm)	762-653	775757626536
	<b>Lid for Fracture Bedpan</b> —	Fracture Bedpan	10" x 8" + tab (25.5 x 20.6 cm +tab)	—	200	19.5" x 9.65" x 10.5" (49.5 x 24.5 x 26.5 cm)	762-654	775757626543
	<b>Large Bowl</b> Versatile vessel, can be used in many settings: large emesis basin, irrigation, etc. Cut-out facilitates proximity to patient. Easy to handle: flat bottom surface and cut-out shape. Lightweight and easy to carry.	—	10.8" x 3.75" (27.5 x 9.5 cm)	2 L	200	26.77" x 21.9" x 11.26" (68 x 55.6 x 28.6 cm)	762-678	775757626789
	<b>Lid for Large Bowl</b> —	Large Bowl	7.8" + tab (19.8 cm + tab)	—	200	19.5" x 9.1" x 8.27" (49.5 x 23 x 21 cm)	762-679	775757626796
	<b>Small Bowl</b> Versatile vessel, can be used in many settings: small emesis basin, irrigation, etc. Cut-out facilitates proximity to patient. Easy to handle: flat bottom surface and cut-out shape. Easy to carry and light.	—	9.1" x3.24" (23 x 8.2 cm)	1.3 L	300	24.2" x 9.25" x 17.5" (61.5 x 23.5 x 44.5 cm)	762-676	775757626765
	<b>General Purpose Bowl</b> Large versatile vessel, can be used as a washing bowl (without detergent), as a commode basin / liner, large emesis basin and many other uses. Easy to handle: flat bottom surface and top lip. Lightweight and easy to carry.	—	12" x4" (30 x 10 cm)	3 L	100	21.9" x 12" x 12.2" (55.5 x 30.5 x 31 cm)	762-700	775757627007
	<b>Kidney Dish – Emesis</b> Standard hospital kidney dish, versatile vessel. Easy to handle: flat bottom and top surfaces with small rim. Easy to carry.	—	9.85" x 5.5" x 1.77" (25 x 14 x 4.5 cm)	700 ml	300	15.75" x 11.6" x 10.6" (40 x 29.5 x 27 cm)	762-650	775757626505
	<b>Female Urinal</b> Standard female urinal. Practical to avoid lifting a patient that is at risk of further injury, onto a bedpan. Comfortable, soft surface, smooth and warm to the touch. Will not adhere to moist skin. Easy to handle. Will stay in place. Lightweight and easy to carry. Can be weighed for volume measurement.	—	9.85" x 4.33" x 4.25" (25 x 11 x 10.8 cm)	900 ml	120	34.65" x 29.5" x 17.9" (88 x 75 x 45.5 cm)	762-655	775757626550
	<b>Male Urinal</b> Designed for patient comfort and caregiver ease of use. Comfortable, soft surface, smooth and warm to the touch. Will not adhere to moist skin. Easy to handle: flat bottom and end surfaces. Will stay in place. Lightweight and easy to carry. Can be weighed for volume measurement.	—	7.87" x 4.53" x 5" (20 x 11.5 x 12.8 cm)	900 ml	110	26.2" x 16" x 27.6" (66.5 x 40.5 x 70 cm)	762-660	775757626604
	<b>Urine Specimen</b> Safe and practical to collect liquid without the fear of being in contact with liquids or of contamination. Can be used with a pail, in a commode chair. Comfortable, soft surface, smooth and warm to the touch. Will not adhere to moist skin. Easy to handle: flat bottom surface and large lip. Graduated to 0.7 L. Lightweight and easy to carry. Can also be weighed for a volume measurement.	—	11.9" x 8.2" x 3.75" (30.27 x 20.7 x 9.5 cm)	1 L	200	17.7" x 26.4" x 17.3" (45 x 67 x 44 cm)	762-701	775757627014