

NEW!

AquaSense®

For a safer bathroom

Raised Toilet Seat



1 Maximum security and support



2 Easy lateral access and secure grip



3 Maximum space and easy access

3 Configurations in **1**

MAX.
300 lb
(136 kg)

Weight capacity



Easy to clean



Adds 4" of height

NEW!


AquaSense®

For a safer bathroom


Raised Toilet Seat

Makes it easier to get on and off the toilet by adapting to the user's specific needs.


- 3 easy-to-assemble configurations for maximum versatility.
- Molded seat for greater comfort.
- Raises toilet seat height by 4" (10 cm) for easier access and use.
- Attaches securely to most standard and elongated toilets for maximum usability.
- Large opening; 10" x 8½" (25.4 cm x 21.6 cm) for optimal fluid flow and less mess.
- Easy-to-clean polyethylene for easy clean ups.
- Corrosion resistant for increased durability.
- Comfortable, easy-grip foam padded armrests.
- Easy-to-turn adjustment knob.
- 300 lb weight capacity for greater stability.
- Limited lifetime warranty.
- Parts available, because AMG is committed to quality and great customer service.



Use both armrests for more security and support getting up and down.



Use one armrest for easy lateral access and a secure grip.



Use without armrests for maximum space and easy access.



Corrosion resistant



Easy to clean



Adds 4" of height



Weight capacity



6.2 lb (2.8 kg)

Product Weight

20.5" - 17" (52 - 43 cm)

Seat Width

16" (40.6 cm)

Seat Depth

18" - 14" (45.7 - 35.6 cm)


Width Between Arms

4" (10 cm)

Raised Height



AMG No.	Description	Retail Box Dimensions (W x D x H)	Unit UPC
770-618	3-in-1 Raised Toilet Seat	18.1" x 7.9" x 21.2" (46 cm x 20 cm x 54 cm)	7 757571 70618 4



The new AquaSense® 3-in-1 Raised Toilet Seat comes in a versatile retail box that features multiple display options.

To learn more about this AquaSense® Raised Toilet Seat or any other AquaSense product, contact your AMG representative, call 1-800-363-2381 or visit: www.myAquaSense.com

AMG AMG Medical Inc.
Montreal, QC H4T 1V5
www.amgmedical.com

REV 1 0209