

MedPro®

The Professional Choice



Balance IV Stand

More strength, more stability in less space.




A stronger, more stable pole

- Wider base for improved stability*
- Lower center of gravity
- Deeper tipping point – tested with up to 30 lbs of added weight (both loaded and unloaded)[†]
- 3" casters for easier glide on different surfaces

Efficient design

- Up to 24% lighter than similar poles[†]
- Less packaging: more environmentally conscious
- Pole available in 2 finishes : Chrome and Stainless Steel
- Optional patient handle and drainage bag hook



Ordering Guide

AMG n°	Description	Pole Finish
775-830	MedPro® Balance IV Pole	Chrome
775-840	MedPro® Balance IV Pole	Stainless Steel
B03-220	Optional Patient Handle	Chrome
B03-230	Optional Drainage Bag Hooks	Chrome

Product Specifications

Caster size:	3" (7.5 cm)	Base diameter:	24" (61 cm)
Height range:	68" to 95" (168 cm to 241 cm)	Total weight:	20.3 lbs (9.2 kg)

*Compared to AMG IV Pole and competitors [†]Internal data

For more information about MedPro® products, contact your AMG representative or our customer service at 1-800-363-2381.

AMG AMG Medical inc.
www.amgmedical.com

REV 2 0310