

MedPro®

The Professional Choice



Aluminum Crutches

- Lighter than wood, they're manufactured from lightweight anodized aluminum to provide professional quality comfort, performance and durability.
- Double-extruded center tube for additional weight-bearing strength up to 300 lbs (136 kg).
- A silencing collar reduces noise and vibration while high impact ABS plastic handles and armrests provide unparalleled security.
- Dual brass push-buttons facilitate adjustments while approximate patient height settings are labeled on the center tube to within 1" increments.
- A curved stair deflector prevents the crutch from becoming lodged and increases stress tolerance at all ground contact points.
- Pre-packaged with accessories made from dry natural rubber which contain denatured proteins to greatly reduce the risk of an allergic reaction.
- MedPro® crutches are backed by a one year limited warranty and our dedicated customer service representatives are always ready to take your call.

Underarm Pads

Thick cushioning, internal ribbing and reinforced ends reduce pressure and underarm fatigue.

Hand-Grips

Washable hand grips feature natural sponge cushioning to absorb impact shock and are perspiration resistant.

Crutch Tips

Ribbed contoured tips provide exceptional traction while a built in metal ring protects against premature wear.



Ordering Guide

AMG n°	Crutch Length	Patient Height	Format	UPC (case)
720-781	52" to 60"	5' 10" to 6' 6"	Pair	7 7575720781 0
720-786	44" to 52"	5' 2" to 5' 10"	Pair	7 7575720786 5
720-791	36" to 45"	4' 6" to 5' 2"	Pair	7 7575720791 9
720-793	31" to 37"	4' to 4' 6"	Pair	7 7575720793 3
720-796	27" to 32"	3' 7" to 4'	Pair	7 7575720796 4

AMG n°	Crutch Accessories	Format	UPC (case)
725-100	Tips, Medium	Pair	7 7575725100 4
725-110	Tips, Large	Pair	7 7575725110 3
725-120	Tips, Extra-large	Pair	7 7575725120 2
725-140	Hand Grips, Open	Pair	7 7575725140 0
725-150	Hand Grips, Closed	Pair	7 7575725150 9
725-160	Underarm Pads	Pair	7 7575725160 8

An online fitting guide can be found on our website at www.amgmedical.com

For more information about MedPro® products, contact your AMG representative or call 1-800-363-2381.



AMG Medical Inc
Montreal, QC H4T 1V5
www.amgmedical.com

REV 1 0109