

Physio Logic - Product Comparison Sheet



BPs

EssentiaA
106-930

SensA
106-910

LuminaA
106-934

SignA
106-937

	EssentiaA 106-930	SensA 106-910	LuminaA 106-934	SignA 106-937
Display Screen	84 mm x 55 mm	76 mm x 57 mm Backlit	102 mm x 69 mm Backlit	105 mm x 55 mm
Memory	Multi-User (4 x 30)	Dual-User (2 x 90)	Dual-User (2 x 60)	Dual-User (2 x 60)
Hypertension Risk Indicator	Yes	Yes	Yes	Yes
Irregular Heartbeat Detection	No	Yes	Yes	Yes
Preset Reminders	No	No	No	Yes
Time & Date	Yes	Yes	Yes	Yes
Room Temperature	No	No	No	Yes
AC Power	With Optional AC Adaptor	With USB	With Optional AC Adaptor	With Optional AC Adaptor
Universal Cuff	Yes (22 - 42 cm)	Yes (22 - 42 cm)	Yes (22 - 42 cm)	Yes (23 - 43 cm)
Warranty	2 Years	5 Years	2 Years	2 Years