


700-869 (MIDNIGHT-BLUE) HUGO TRANSPORT CHAIR 19" WIDE



**700-869 (MIDNIGHT-BLUE)
HUGO TRANSPORT CHAIR 19" WIDE**

Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	C48-00076	Blue Folding Handle & Quick Fold Mechanism - L & R	1	1	PR
2.	C48-00077	Handle Release Lever, Black - 2 ea / set	1	1	ST
3.	C48-00078	Full Black Nylon Arm Pad & Hardware – 2 ea / set	1	1	ST
4.	C48-00061	Seat Strap - 2 ea /set	1	1	ST
5.	C48-00079	Black Nylon Seat Upholstery 19"x15.5" & Hardware	1	1	EA
6.	C48-00081	Blue Swing-Away Footrest - L & R	1	1	PR
7.	C48-00082	Black Footrest Release Mechanism - L & R	1	1	PR
8.	C48-00066	Safety Belt, Black	1	1	EA
9.	C48-00084	Hugo Black Nylon Back Upholstery 19"x16" & Hardware	1	1	EA
10.	C48-00069	2 Front & 2 Rear Black Seat Guides	1	1	ST
11.	C48-00070	Wheel Lock Assembly - L & R	1	1	PR
12.	C48-00087	2 Front + 2 Rear End Plugs	1	1	ST
13.	C48-00085	8" Black Rear 6 Spokes Wheel & Hardware	2	1	EA
14.	C48-00086	Set of Knobs & Hardware – 6 ea / set	1	1	ST
15.	C48-00088	Front Fork & Hardware	2	1	EA
16.	C48-00089	8" Black Front 6 Spokes Wheel & Hardware	2	1	EA