


770-507 AQUASENSE ADJUSTABLE BATH BENCH - YELLOW



Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	A11-403	Seat, Blow-Molded Plastic, Yellow (505x305x35mm)	1	1	EA
2.	B11-823	Hardware Set, Yellow (5 Triangular Knobs M6x34mm)	1	1	ST
3.	B11-811	Push Button dia:8mm & Spring - 4ea / set	1	1	ST
4.	A11-213	Anti-Rattle for dia:28.6mm, Yellow - 4 ea / set	1	1	ST
5.	B11-013	Extension Leg, Yellow (Aluminum Tube dia:28.6x250mm, Anti-Rattle & Rubber Tip) - 4 ea / set	1	1	ST
6.	A11-483	Rubber Tip for dia:28.6mm, Yellow - 4 ea / set	1	1	ST
Not Shown	B11-851	Hardware Set for Bath Bench w/out Back (5 Cross-Cup Screws M6x34mm, 5 Curved Washers) - Accessory	1	1	ST